

## MICROSCOPE OBJECTIVE

MLOptic not only provide standard microscope objectives to fit your existing system, but also provide high-end customized objectives to help your system achieve a higher level of performance.

### ■ M Plan Apo Series

#### M PLAN APO

HIGH N.A.  
FLAT FIELD / APO  
VIS

1

2

3

4

5

6

7

#### M PLAN APO L

LONG WORKING DISTANCE  
HIGH N.A.  
FLAT FIELD / APO  
VIS

#### M PLAN APO HR

HIGH RESOLUTION  
FLAT FIELD / APO  
VIS

#### M PLAN APO NIR

HIGH RESOLUTION  
FLAT FIELD / APO  
VIS-NIR

#### M PLAN APO NIR L

LONG WORKING DISTANCE  
HIGH N.A.  
FLAT FIELD / APO  
VIS-NIR


#### M PLAN APO NIR HR

LONG WORKING DISTANCE  
HIGH RESOLUTION  
FLAT FIELD / APO  
VIS-NIR


#### M PLAN APO NUV

LONG WORKING DISTANCE  
HIGH N.A.  
FLAT FIELD / APO  
VIS-NUV






## MICROSCOPE OBJECTIVES – STANDARDIZED SERIES

	Mag.	f' (mm)	NA	WD (mm)	FOV (mm)	DOF (μm)	Mount
	2X	100	0.055	≥34	φ5.5	91	M26X0.706
	5X	40	0.14	≥34	φ2.23	14	M26X0.706
	10X	20	0.28	≥33.5	φ1.1	3.5	M26X0.706
	20X	10	0.42	≥20	φ0.55	1.6	M26X0.706
	50X	4	0.55	≥13	φ0.22	0.9	M26X0.706

## MICROSCOPE OBJECTIVES – CUSTOMIZED SERIES

#	Mag.	f' (mm)	NA	WD (mm)	FOV (mm)	DOF (μm)	PD (mm)	Mount	
1	20X	10	0.75	≥1	∅1.25	0.98	60	M26X0.706	
2	20X	10	0.8	≥0.8	∅1.57	0.86	60	M26X0.706	
3	40X	5	0.95	≥0.21	∅0.625	0.61	60	M26X0.706	

■ **MLOPTIC Corp. Headquarters**  
 ADD: 398 Pugang St., Jiangning Development Zone Nanjing, 211102 China  
 TEL: 0086-25-84436288 FAX: 0086-25-84433188 52728163 E-Mail: [info@mloptic.com](mailto:info@mloptic.com)  
 WEBSITE: [www.mloptic.com](http://www.mloptic.com)

■ **MLOPTIC (THAILAND) CO.,LTD.**  
 Contact Name: Mike Young  
 ADD: No. 700/154 Moo 1 Amata Nakorn Industrial Estate ,Ban Kao Subdistrict, Phan Tong District,  
 Chonburi Province Thailand  
 Mobile: +66 0910 5841 21  
 E-Mail: [mike@mlthailand.com](mailto:mike@mlthailand.com)

